

Technical data



Oticon Cara™ BTE T 105

Oticon Cara™ BTE T 105 comes with an undamped hook and Corda miniFit. It features a disposable battery. It is built on the Sirius™ platform, powered by Oticon BrainHearing™ technology and AI-driven sound processing. It supports

Bluetooth® LE Audio, Bluetooth Low Energy, Auracast™ broadcast and Fast Pair. It enables hands-free communication and direct streaming for iPhone, iPad, Vision Pro, Apple Watch, Mac and select Android™ devices.

Hook



BTE T 105

Corda miniFit 1.3 mm



BTE T 105

Corda miniFit 0.9 mm



BTE T 105

Technical Features

- › Auracast™ broadcast
- › Bluetooth® LE Audio
- › Bluetooth® Low Energy
- › Fast Pair
- › Telecoil
- › NFMI (Near-Field Magnetic Induction)
- › Double push-button
- › IP68 rated

Accessories

- › Oticon Companion app
- › ConnectClip
- › TV Adapter 3.0
- › EduMic
- › Phone Adapter 2.0
- › Hook and Corda miniFit (thin tube)
- › Size 13 zinc air battery

For information on compatibility, please visit www.oticon.co.uk/products-nhs/compatibility-guide.

Operating conditions

Temperature: +1°C to +40°C (+34°F to +104°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

Transportation and storage conditions

Temperature and humidity shall not exceed the mentioned limits during transportation and storage.

Transportation

Temperature: -25°C to +60°C (-13°F to +140°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

Storage


Temperature: -25°C to +60°C (-13°F to +140°F)
Humidity: 5% to 93% relative humidity, non-condensing
Atmospheric pressure: 700 hPa to 1060 hPa

WARNING: No modification of this equipment is allowed.

Apple, the Apple logo, iPhone, iPad, Vision Pro, Apple Watch, Mac and the Mac logo are trademarks of Apple Inc., registered in the U.S. and other countries. Use of the Made for Apple badge means that an accessory has been designed to connect specifically to the Apple product(s) identified in the badge and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Google, Android, and related marks and logos are trademarks of Google LLC. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by the Bluetooth SIG. Any use of such marks by Demant is under license. Other trademarks and trade names are those of their respective owners.

Fitting ranges

Oticon Cara



Feature overview

	Oticon Cara
Speech understanding & listening ease	
MoreSound Intelligence™ 3.0	Level 2
Environment classifier	5 configurations
Virtual Outer Ear	2 configurations
Spatial Balancer	60%
Neural Noise Suppression, Difficult / Easy	10 dB / 4 dB
Sound Enhancer	2 configurations
Wind & Handling Stabilizer	•
MoreSound Amplifier™ 3.0	•
SuddenSound Stabilizer	5 configurations
MoreSound Optimizer™	•
Feedback shield	•
Spatial Sound™	4 estimators
Soft Speech Booster	•
Frequency lowering, Speech Rescue™	•
Sound quality	
Clear Dynamics	•
Better-Ear Priority	•
Fitting Bandwidth ¹	8 kHz
Power Bass (streaming)	•
Processing Channels	48
Personalisation & optimised fitting	
Fitting Bands	20
Multiple Directionality options	•
Adaptation Management	•
Fitting Formulas	VAC+, NAL-NL1/NAL-NL2, DSL v5
Connecting to the world	
Oticon Companion app	•
Bluetooth® LE Audio ²	•
Auracast™ broadcast ²	•
Hands-free communication ²	•
Direct streaming ²	•
ConnectClip	•
EduMic	•
Remote Control 3.0	•
TV Adapter 3.0	•
Phone Adapter 2.0	•
Tinnitus SoundSupport™	•
CROS/BiCROS support	•

1) Bandwidth accessible for gain adjustments during fitting

2) Available on select devices. For more information, please visit www.oticon.co.uk/products-nhs/compatibility-guide

Oticon Cara BTE T 105

Ear Simulator



Measured according to IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2022, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010



Technical information
Omnidirectional mode is used unless otherwise stated.

105 hook (undamped) / Corda miniFit 1.3
Acoustic input: 60 dB SPL ———
Magnetic input: 31.6 mA/m - - - - -

Corda miniFit 0.9
Acoustic input: 60 dB SPL ———
Magnetic input: 31.6 mA/m - - - - -



OSPL90, Peak (dB SPL)	139	134	130
OSPL90, 1600 Hz (dB SPL)	130	123	116
OSPL90, HFA (dB SPL)	133	128	121
Full-on gain, Peak (dB) ¹	72	70	67
Full-on gain, 1600 Hz (dB) ¹	63	56	49
Full-on gain, HFA (dB) ¹	67	62	55
Reference test gain (dB)	56	48	42
Frequency range (Hz)	<100-7500	<100-7500	<100-7500
Telecoil output, 1 mA/m field (1600 Hz) (dB SPL)	94	87	81
Telecoil output, 10 mA/m field (1600 Hz) (dB SPL)	114	107	101
Total harmonic distortion (Input 70 dB SPL), 500 Hz (%)	7	<2	<2
Total harmonic distortion (Input 70 dB SPL), 800 Hz (%)	5	<2	<2
Total harmonic distortion (Input 70 dB SPL), 1600 Hz (%)	<2	<2	<2
Equivalent input noise level, Omni (dB SPL)	19	23	28
Battery consumption, Typical (mA) ²	2.1	2.1	2.1
Battery consumption, Quiescent (mA) ²	2.0	2.0	2.0
Battery life, artificial measurement, hours ³	150	150	150
Expected battery life, hours (battery size 13 - IEC PR4B) ⁴	85-105	85-105	85-105

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response but without influence of feedback.

2) Battery current is measured after a settling time of minimum 3 minutes.

3) Based on the standardised battery consumption measurement. The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.

4) Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, incl. direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).

Warning to the hearing aid dispenser

The maximum output capability of the hearing aid may exceed 132 dB SPL (2cc coupler). Special care should be exercised in selecting and fitting the hearing aid, as there may be risk of impairing the remaining hearing of the hearing aid user.

Oticon Cara BTE T 105

2CC Coupler



Measured according to ANSI S3.22-2024,
IEC 60118-0:2022 and IEC 60318-5:2006



Technical information
Omnidirectional mode is used
unless otherwise stated.

105 hook (undamped) / Corda miniFit 1.3
Acoustic input: 60 dB SPL ———
Magnetic input: 31.6 mA/m - - - - -

Corda miniFit 0.9
Acoustic input: 60 dB SPL ———
Magnetic input: 31.6 mA/m - - - - -



OSPL90, Peak (dB SPL)	132	130	126
OSPL90, HFA (dB SPL)	126	119	113
Full-on gain, Peak (dB) ¹	66	66	63
Full-on gain, HFA (dB) ¹	61	54	47
Reference test gain (dB)	50	43	36
Frequency range (Hz)	<100-7000	<100-7300	<100-7300
Telecoil output, HFA SPLITS L/R (dB SPL)	109	100	94
Telecoil output, Full-on HFA-SPLIV (dB SPL)	110	105	98
Total harmonic distortion (Input 70 dB SPL), 500 Hz (%)	4	<2	<2
Total harmonic distortion (Input 70 dB SPL), 800 Hz (%)	3	<2	<2
Total harmonic distortion (Input 65 dB SPL), 1600 Hz (%)	<2	<2	<2
Total harmonic distortion (Input 60 dB SPL), 3200 Hz (%)	<2	<2	<2
Equivalent input noise level, Omni (dB SPL)	15	20	26
Battery consumption, Typical (mA) ²	2.4	2.5	2.4
Battery consumption, Quiescent (mA) ²	2.0	2.0	2.0
Battery life, artificial measurement, hours ³	130	125	130
Latency, (ms)	8.4	8.4	8.4
Expected battery life, hours (battery size 13 - IEC PR48) ⁴	85-105	85-105	85-105

1) Measured with the gain control of the hearing aids set to their full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response but without influence of feedback.

2) Battery current is measured after a settling time of minimum 3 minutes.

3) Based on the standardised battery consumption measurement. The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.

4) Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, incl. direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).

Warning to the hearing aid dispenser

The maximum output capability of the hearing aid may exceed 132 dB SPL (2cc coupler). Special care should be exercised in selecting and fitting the hearing aid, as there may be risk of impairing the remaining hearing of the hearing aid user.

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



SBO Hearing A/S
Kongebakken 9
DK-2765 Smørum
Denmark

Headquarters
Oticon A/S
Kongebakken 9
DK-2765 Smørum
Denmark