



OTICON | Play SI

Fitting range*

miniRITE R DSL Fitting ranges



miniBTE R DSL Fitting ranges



*Fitting range is based on Play SI 1. Fitting range for Play SI 2 is limited to 8kHz. Details available in Technical Data sheets.

**Mold includes MicroShell Detect, MicroMold & LiteTip.

Oticon Play SI™

miniRITE R and miniBTE R




Oticon Play SI miniBTE R offers both fitting level 85 and 105 giving you more flexibility during fittings and your patients more amplification when needed. This style also features a size 13 lithium-ion battery and telecoil.

Oticon Play SI miniRITE R is modern and sleek. It features both a telecoil and a size 13 lithium-ion battery that gives hearing aid users more power between charges.

miniFit Detect receivers

Our self-calibrating innovation gives precise gain to the user and valuable information to you about fitting, thanks to an embedded memory. Select between 4 receivers with different output performance. The 105 receiver is only available as a custom MicroShell Detect.




Hook and Thin tube for miniBTE

Select between miniBTE hook with our custom BTE molds or Corda miniFit thin tubes with one of our miniFit earpieces.



Earpieces

miniFit Detect earpieces include standard domes and custom molds. LiteTip is available in acrylic or OtoTherm™ material. MicroMold is available in acrylic, soft silicone, or OtoTherm material. MicroShell Detect is only available in acrylic.



IP68

Proven IP68 water and dust resistant hearing aids.



	Oticon Play SI 1	Oticon Play SI 2
Speech understanding & listening ease	MoreSound Intelligence™ 3.0	Level 1
	Sensor technology	✓
	Environment classifier	5 configurations
	Virtual Outer Ear	3 configurations
	Spatial Balancer	100%
	Neural Noise Suppression Difficult/Easy	12 dB / 6 dB
	Sound Enhancer	3 configurations
	Wind & Handling Stabilizer	✓
	MoreSound Amplifier™ 3.0	✓
	SuddenSound Stabilizer	6 configurations
	MoreSound Optimizer™	✓
	Feedback shield	✓
	Spatial Sound™	4 estimators
	Soft Speech Booster	✓
Sound quality	Frequency lowering, Speech Rescue™	✓
	Clear Dynamics	✓
	Better-Ear Priority	✓
	Fitting Bandwidth ¹	10 kHz
	Power bass (streaming)	✓
Personalization & optimizing fitting	Processing channels	64
	Fitting bands	24
	REM Autofit ²	Verifit® LINK, IMC 2
	Pediatric Fitting Mode	✓
	DSL Fitting Range ³	✓
Designed for children	Fitting formulas	DSL v5, NAL-NL 1/ NAL-NL 2, VAC+
	EduMic	✓
	LED	✓
	Biologically safe	✓
	Hearing aid color options	12
	Hook color options	6
	Bluetooth® LE Audio ⁴	✓
	Auracast™ broadcast ⁴	✓
	Hands-free communication ⁴	✓
	Fast Pair	✓
	Oticon Companion app	✓
	Tinnitus SoundSupport™	✓
	CROS/BiCROS support	✓

¹ Bandwidth accessible for gain adjustments during fitting.² IMC 2: Inter Module Communication 2³ Available in Technical Data and marketing materials⁴ For more information on compatible phones, please visit oticon.com/support/compatibility-guide.

OTICON | Play SI

Product comparison



oticon
life-changing technology

Apple, the Apple logo, iPhone®, iPad®, Mac®, and the Mac logo are trademarks of Apple Inc., registered in the U.S. and other countries. Use of the Made for Apple badge means that an accessory has been designed to connect specifically to the Apple product(s) identified in the badge and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Android™ is a trademark of Google LLC. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by Bluetooth SIG. Any use of such marks by Demant is under license. Other trademarks and trade names are those of their respective owners.

Oticon and Play SI are trademarks of Oticon A/S, registered in the United States and other jurisdictions. All rights reserved.

