

E) Available CROS/BiCROS and Bimodal solutions

CROS solution



microRIE CROS Transmitter

BTE CROS Transmitter
(not available for ReSound Vivia)

Ensure a high quality and open sound picture for users with single-sided deafness

Cochlear™ bimodal solution




Smart Hearing Alliance

If a hearing aid and a hearing implant are needed, ReSound and Cochlear's solutions work seamlessly together¹

F) Choose how you want to control and connect

Supports all ReSound hearing aids except CIC styles

Control at your fingertips




ReSound Smart 3D™ app

Access to controls at all times and Auracast broadcast audio with the world's first Auracast Assistant for hearing aids.



ReSound Remote Control

Control at your fingertips



Crystal clear phone calls

Hands-free calls




Bluetooth® Low Energy Audio

Connecting hands-free to Bluetooth 5.3 compatible Android devices²

Hands-free calls for iPhone and iPad

Connecting hands-free to iPhone and iPad³

Streamed calls



Direct audio streaming

From iPhone, iPad and compatible Android™ devices⁴

Not supported with ReSound Key™ tech level 2

Bring harmony to TV time



ReSound TV Streamer 2

For hearing great TV audio



ReSound TV-Streamer+

For immediate TV speech clarity at your preferred volume and with Auracast broadcast support

Compatible with ReSound Vivia, ReSound Savi™ and ReSound Nexia

Hear conversations more clearly



ReSound Multi-Mic+

Cut through the noise now with the world's first remote mic with Auracast™

Compatible with ReSound Vivia, ReSound Savi and ReSound Nexia



Telecoil

Available in Multi-Mic+, RIE 13 and BTE styles



Discover more on gnhearing.com

Worldwide Headquarters GN Hearing A/S DK-2750 Ballerup Denmark Tel.: +45 4575 1111 gnhearing.com CVR no. 55082715	United Kingdom GN Hearing UK Ltd. Unit 13 Talisman Business Centre Bicester Oxon OX26 6HR Tel: +44 1869 352800 gnhearing.com	Australia GN Hearing Australia Pty Ltd Suite 1001, Level 1 1 Julius Avenue North Ryde, NSW 2113 Australia Tel: (free) 1800 658 955 gnhearing.com	New Zealand GN Hearing New Zealand Limited Ground Floor, North Entrance 4 Fred Thomas Drive Takapuna, Auckland, 0622 Tel.: (free) 0800 900 126 gnhearing.com
--	--	--	---

1 ReSound Vivia and ReSound Savi are not bimodally compatible. Cochlear and GN continue to work together to bring state of the art bimodal solutions to Hearing Care Professionals and users under the Smart Hearing Alliance

2 Compatible with ReSound Vivia, ReSound Savi and ReSound Nexia

3 Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation), iPad (10th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later.

4 See compatibility on gnhearing.com/products/compatibility

M2017036B-25-03-Rev.A

© 2025 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone and iPad are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast word mark and logos are trademarks owned by Bluetooth SIG, Inc. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. Cochlear, Hear now, And always, Nucleus, Kanso, Baha, the elliptical logo, and marks bearing an ® or ™ symbol, are either trademarks or registered trademarks of Cochlear Limited or Cochlear Bone Anchored Solutions AB (unless otherwise noted).


For professionals

ReSound GN


GN Bringing People Closer

... starts with understanding your daily hearing challenges, surroundings and preferences. Because everyone spends their time differently, understanding your everyday is critical when choosing a solution that will support you in your favourite activities.


1. How important is it for you to improve your hearing right now?




2. How does your hearing limit you in doing what you want to do?




3. How does your hearing ability affect your social life?



4. How is your hearing holding you back in your professional life?



The figure depicts the sounds of daily life and where they typically fall relative to pitch and loudness. Frequency (pitch) is measured on the x-axis, while sound intensity (loudness) is measured on the y-axis. Letters on the graph denote common sounds of speech. The banana shaped area shows the pitch and loudness of these common speech sounds.



Very noisy

Noisy

Moderately noisy

Quiet

Essential

Advanced

Premium

The amount of noise you experience in your daily life is one of the deciding factors when choosing a hearing aid. This diagram illustrates how you can benefit from having more advanced features that boost clarity and unblur conversations in noisier environments.

As you move through the tech levels, from Essential to Advanced and Premium, the greater the degree that noise is reduced, speech is highlighted, and sound is crystal clear.

At the ultimate level – Premium – you will experience a vivid sound picture, where you still have access to important noises, but speech is contrasted, enabling you to enjoy greater clarity in your listening experiences.

	Characteristics	Styles		
Receiver-In-Ear (RIE)	- Mild to severe hearing loss - All ear canal types - Hear more naturally - Discreet	Rechargeable microRIE 60S miniRIE 60 61 Premium charger Standard charger Desktop charger		
		Non-rechargeable 61 312 battery 62 13 battery		
Behind-The-Ear (BTE)	- Mild to severe hearing loss - All ear canal types - Easy to use and robust	Rechargeable 77 Desktop charger		
		Non-rechargeable 67 312 battery 77 13 battery		
Power Behind-The-Ear (BTE)	- Severe to profound hearing loss - All ear canal types - Easy to use and robust	Rechargeable 88 Plastic hook 88 Metal hook Desktop charger	Non-rechargeable 88 Plastic hook 88 Metal hook	
Super Power Behind-The-Ear (BTE)		Rechargeable 98 Metal hook 98 Plastic hook Desktop charger	Non-rechargeable 98 675 battery	
Customs In-The-Ear (ITE)/In-The-Canal (ITC)	- Mild to severe hearing loss - Narrow, curvy and standard ear canals - Earbud experience and hear more naturally	Rechargeable ITC ITE Custom charger Devices will be built as small as possible based on your ear anatomy.	Non-rechargeable ITC 312 battery	
Customs Completely-In-Canal (CIC)		Non-rechargeable NW-CIC 10A battery		