For professionals

A portfolio of hearing solutions **tailored to your clients**



Beltone **Portfolio**

Beltone

Welcome more clients with the most intelligent Beltone Portfolio yet

Organic Hearing™ philosophy

Our Organic Hearing philosophy starts with real people and real listening situations. It is a human-centred approach to developing fit-for-purpose in-clinic solutions and innovative technology.

By using the human brain as our reference for how sound is naturally processed, we can deliver hearing aids that sound natural, feel natural and help your clients stay connected.



Sounds natural

- Top-rated for hearing in noise[1]
- 88% better accuracy in real-time classification and feature engagement[2]
- 80% of participants experience better audibility of the intended talker[3]



Feels natural

- Comfortable all-day wear
- All-weatherproof designs
- All-day battery life



Connects naturally

- Made for Bluetooth® Low Energy (LE) Audio and Auracast™
- Clear hands-free calls for iPhone, iPad and compatible Android™ devices*
- Wireless accessories to stream high-quality audio and speech

Our human approach to AI

We believe AI should empower people, not replace them. That belief sits at the heart of our unique approach to AI. Inspired by our Organic Hearing philosophy, we use AI to support the brain's natural way of processing sound, delivering the best possible listening experience for each and every Beltone hearing aid user.

1 Schumacher & Groth (2025), 2 Compared to closest Premium brand. Groth & Cui (2026), 3 Compared to other Premium hearing aid brands with DNN denoising. Jespersen & Rubachandran (2026).

* Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation), iPad (10th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later. Bluetooth LE Audio works with Beltone Illuminate, Beltone Envision, Beltone Boost Max S, Beltone Commence, Beltone Serene and compatible mobile devices with Bluetooth 5.3 or later. See compatibility on gnhearing.com/products/compatibility

Hearing solutions to delight your clients

Premium — **Beltone Illuminate**

RIE and Custom styles New
Mild to profound hearing loss



CIC Non-wireless · ITC Wireless · ITC RHI · ITE RHI · microRIE 62S RHI · RIE 63 312 battery · RIE 64 13 battery

Premium — **Beltone Envision**

RIE styles
Mild to profound hearing loss



microRIE 62S RHI · RIE 63 312 battery · RIE 64 13 battery

Premium — **Beltone Boost Max S**

BTE styles
Severe to profound hearing loss



Super Power BTE 95
Metal hook / Plastic hook
RHI

Premium — **Beltone Serene**

RIE, BTE and Custom styles
Mild to profound hearing loss



CIC Non-wireless · ITC Wireless · ITC RHI · ITE RHI · microRIE 62S RHI · RIE 63 312 battery · RIE 64 13 battery · BTE 72 RHI · BTE 76 RHI · Power BTE 86 RHI

Essential — **Beltone Commence Beltone Boost SP**

RIE and Custom styles
Mild to profound hearing loss



CIC Non-wireless · ITC Wireless · miniRIE 62 RHI · RIE 63 312 battery · RIE 64 13 battery



Super Power BTE 95*

Receivers, thin tubes and domes overview

Auracast™ and streaming

Connectivity and accessories

Feature overview

RHI = Rechargeable Hearing Instrument

* Beltone Boost SP

Everything you need **in one place**

New

Beltone Illuminate™

Bring hearing clarity to life even in complex listening environments with our new Premium hearing aid family. Featuring our new AI platform and deep neural network (DNN) technology, Beltone Illuminate delivers exceptional clarity to your clients' auditory world.



Beltone Commence™

Our Essential family with RIEs and Customs gives your clients everything they need, and a little more. Offer your clients best-in-class hearing all around, hearing comfort and Bluetooth connectivity.



Beltone Boost™ Max S

Introducing super-small, super-powerful rechargeable hearing aids with all-day battery for your clients with severe to profound hearing loss. Boost hearing in noise, amplification, feedback management, connectivity, and more with our latest hearing innovations.



Bring your best care to every client


Introducing the fully upgraded Beltone Portfolio, thoughtfully designed for your business, your practice and your clients.



For your business

Exceptional hearing care goes beyond innovative technologies. The Beltone Portfolio delivers life-changing clarity and the support you need to keep making a difference every day.

- **More options for your clients**
  With our Premium and Essential hearing aid families, your clients can immediately get the latest hearing technology in their ideal style without having to compromise.
- **Specialist solutions in all tech levels**
  Bring the latest hearing technology and styles even to clients with profound hearing loss or single-sided deafness with our leading Super Power BTE and CROS offerings.
- **An upgrade on upselling**
  Help your clients hear the value of different hearing aid options with our flexible demo devices covering all tech levels and our latest families.



For your practice

With Beltone Solus Max™, we continue to support your preferred way of fitting and equip you with the tools to help reduce returns and give you more control over inventory and service.

- **Easier, faster, better fittings**
  Beltone Solus Max 2.4 introduces a simplified navigation with an intuitive layout, making it easier to find features. It is also the first fitting software to feature NAL-NL3™.
- **Time-saving solutions for more care**
  Beltone hearing aids come with unified hearing aid parts, including aligned receivers and AllFit domes that are fully backward compatible across all families.
- **Achieve up to 40% fewer service returns**[1]
  Reduce interruptions with on-the-spot exchangeable mic filter replacements and keep the hearing aids on your clients' ears where they belong.



For your clients

Bring the proven benefits of our AI-powered hearing technology to more clients.

- **Life sounds natural automatically**
  Offer brand-new hearing solutions featuring our new AI platform across more programs and tech levels, so your clients can enjoy a more natural hearing experience with less listening effort.
- **Bring best-in-class clarity to more**
  Our leading AcoustIQ™ sound-sensing technology and our unique beamforming strategies are available across all tech levels and price points, bringing advanced sound performance in the smallest styles to more of your clients.
- **All-inclusive connectivity with Auracast**
  Beltone wireless hearing aids enable your clients to connect with Auracast broadcast audio and Bluetooth® LE Audio for the clearest streamed experiences* in larger venues and in smaller conversations.

1 GN Proprietary data on file (2025)

* Bluetooth SIG (2025)

New Beltone Illuminate

Make the world brighter for more people

Enable your clients to experience sound in a whole new light with Beltone Illuminate. Our new Premium family introduces True Directionality*, a new automatic feature comprising 3 types of hearing technology that work in harmony to deliver exceptional clarity in noise.



AcoustIQ

Responds to your surroundings instantly

4-mic Beamforming

Follows your focus for the clearest speech

DNN Denoising

Activates intelligent denoising when needed

VALIDATED

Hörzentrum Oldenburg confirms:

Beltone Illuminate delivers the best audibility advantage compared to other DNN denoising solutions[1] and low listening effort in noise.[2]

1 Jespersen, Heeren & Quilter (2026)
2 Jespersen, Heeren & Quilter (2025)

Exceptional **clarity in noise**

Technology that shines bright

Beltone Illuminate is recognised not only by prestigious institutions, but by the people who use it every day.



Most accurate real-time sound sensing

AcoustIQ delivers real-time sound sensing that adapts accurately to your clients' situation, giving them greater access to the full sound environment.

88% greater accuracy in real-time classification and feature engagement than the closest competitor.[1]



Narrowest 4-mic beamforming

Lets voices shine through by following your clients' focus.

Delivering a **13%** SNR improvement,[2] with **83%** of users achieving superior speech clarity in noise.[3]



Best automatic DNN denoising technology


Reduces sound chaos to bring relief and clarity to your clients' conversations, with fewer distractions.

17x more efficient DNN denoising,[4] resulting in the best audibility of the intended speaker for **80%** of users.[5]



All-day battery you can rely on

Beltone Illuminate hearing aids deliver 24 hours of battery, even with streaming and AI active.[6] The hearing aids are tested in real-life situations, so your clients can be confident they have all the power they need.



Styles that feel natural to wear

Our range of small, all-weatherproof IP68-rated hearing aids offers the all-day comfort users prefer and options to suit different lifestyles – such as the world's smallest AI rechargeable RIE.[6]



Connect with Auracast now

With Beltone Illuminate, your clients can enjoy the clearest streamed sound with Bluetooth® Low Energy (LE) Audio and Auracast** - the universal solution for accessible audio in public and private venues.

1 Groth & Cui (2026)
2 Jespersen & Hashem (2026)
3 Jespersen, Dieu & Rubachandran (2024)
4 Compared to other hearing solutions with AI. Jespersen, Dieu & Rubachandran (2024)
5 Compared to other Premium hearing aid brands with DNN denoising. Jespersen & Rubachandran (2026)
6 GN Proprietary data on file (2026)

* Available in Tech Level 17 for the rechargeable microRIE style
** Bluetooth SIG

New

A brighter **world of sound**

Beltone **Illuminate**

Bring hearing clarity to life with Beltone Illuminate, featuring the most discreet AI hearing aids[1] built on our new AI platform, which adapts automatically with leading precision to illuminate the moments that matter most.

Best audibility advantage in noise. Proven for lower listening effort[2-3]

Tested in real-world scenarios and preferred over competition

Connect with Auracast across the full portfolio*

RIE colours

Sparkling Silver · Sand · Gold · Bronze · Warm Grey · Graphite · Espresso · Deep Black · Navy Blue · Red

Custom colours

Anthracite · Light Beige · Beige · Medium Brown · Brown · Dark Brown

Charging options



Premium charger
Portable with powerbank
Compatible with the microRIE (CFF 8)



Standard charger
Portable
Compatible with the microRIE (CFF 8)



Desktop charger
Easy-to-use
Compatible with microRIE (CFF 8)









Custom charger
Portable for ITC/ITE
(CFF C)

1 GN Proprietary data on file (2026)
2 Jespersen, Heeren & Quilter (2026)
3 Quilter, Heeren & Jespersen (2025))

* Not available for CIC
** All-day battery on one charger is based on an average use time of 20 hours, streaming 50% of the time
*** Only available for ITC
**** Only available for ITE

For mild to profound hearing loss

Style	Features	Fitting range	Max output/gain°
 microRIE 62S-DRWC Rechargeable Available as CROS Transmitter	• Receiver options: MM, LP, MP, HP and UP • All-day battery on one charge** • IP68-rated • Exchangeable mic filters • TapControl with microRIE • Telecoil in 64 • Programmable push button or toggle • CROS: ILU062S-DRWC		**MM**: 115/61 **LP**: 113/56 **MP**: 115/60 **HP**: 120/60 **UP**: 129/71 (Encased)
 RIE 63-DRW 312 battery			
 RIE 64-DRW 13 battery			
 ITC-DWC / ITE-DWC Rechargeable	• Receiver options: LP***, MP, HP and UP • Shell options: ITC, Half-shell**** and Full-shell**** • All-day battery on one charge* • IP68-rated • Programmable push button		**LP**: 114/40*** **MP**: 118/50 **HP**: 121/60 **UP**: 130/70
 ITC-DW 312 battery	• Receiver options: LP, MP, HP and UP • Shell options: ITC, Half-shell and Full-shell • IP68-rated • Exchangeable mic filters • Programmable push button and optional volume control		**LP**: 113/40 **MP**: 115/50 **HP**: 120/60 **UP**: 128/68
 CIC non-wireless 10A battery	• Cabled fitting required • Receiver options: LP, MP and HP • IP68-rated • Optional push button		**LP**: 113/40 **MP**: 115/50 **HP**: 120/60

° ANSI S3.22-2014 / IEC 60118-0:2022 / 2cc coupler

Intelligent hearing
technology and sound

Beltone **Envision™**

Beltone Envision gives beautiful clarity to the noise of life. Help your clients hear with renewed focus and experience sound like never before.

Award winning AI hearing aids for hearing in noise[1-3]

World's smallest AI hearing aid[4]

All-day battery* even when streaming or using True Focus

Colours

Sparkling Silver · Sand · Gold · Bronze · Warm Grey · Graphite · Espresso · Deep Black · Navy Blue · Red

Charging options



Premium charger
Portable with powerbank
Compatible with the microRIE (CFF 8)



Standard charger
Portable
Compatible with the microRIE (CFF 8)






Desktop charger
Easy-to-use
Compatible with the microRIE (CFF 8)

1 Schumacher & Groth (2025)
2 Cui & Groth (2024)
3 Quilter (2021)
4 GN Proprietary data on file

* All-day battery on one charge is based on an average use time of 20 hours, streaming 50% of the time

For mild to profound hearing loss

Style	Features	Fitting range	Max output/gain°
 microRIE 62S-DRWC Rechargeable Available as CROS transmitter	• Receiver options: MM, LP, MP, HP and UP • All-day battery on one charge` • IP68-rated • Exchangeable mic filters • TapControl with microRIE • Telecoil in 64 • Programmable push button or toggle • CROS: ENV162S-DRWC		MM: 115/61 LP: 113/56 MP: 115/60 HP: 120/60 UP: 129/71 (Encased)
 RIE 63-DRW 312 battery			
 RIE 64-DRW 13 battery			

° ANSI S3.22-2014 / IEC 60118-0:2022 / 2cc coupler



Boost your **world of sound**

Beltone **Boost Max S**

Introducing super-small, super-powerful hearing aids for severe to profound hearing loss, featuring our latest innovations to empower those who face the toughest listening challenges.

Made for hearing in noise 2.7dB SNR improvement with Clear Focus[1]

World‘s best feedback management[2]

World‘s smallest rechargeable Super Power hearing aid[3]

Colours

Sparkling Silver · Sand · Gold · Bronze · Warm Grey · Graphite · Espresso · Deep Black · Navy Blue · Red

Charging options





Desktop charger
Easy-to-use
Compatible with the Super Power BTE (CFF 90)

1 For those with severe to profound hearing loss compared to legacy. Jespersen & Skovlund (2025)
2 Compared to leading Super Power solutions. Cui & Groth (2025)
3 Compared to leading rechargeable Super Power solutions. GN Proprietary data on file (2025)

* All-day battery on one charge is based on 28 hours of use – and 20 hours when streaming 50% of the time

For severe to profound hearing loss

Style	Features	Fitting range	Max output/gain°
 Super Power 95-DWHC Metal hook Rechargeable	• All-day battery on one charge* • IP68-rated • TapControl • Programmable push button and volume control • Telecoil • Exchangeable mic filters • Optional plastic hook configurable in FSW		Metal hook: 140/84 Plastic hook: 140/79
 Super Power 95-DWHC Plastic hook Rechargeable			

° ANSI S3.22-2014 / IEC 60118-0:2022 / 2cc coupler



Inspiring a **new generation**

Beltone **Serene™**

Beltone Serene hearing aids connect with Bluetooth LE Audio and Auracast, so your clients can feel more connected, included and social, wherever they are. All while enjoying an incredible sound experience, even in noisy situations, with tiny and comfortable styles.

Industry-best access to surrounding sound[1]

Smallest styles for every need

World's first hearing aids with Bluetooth LE Audio & Auracast broadcast audio[2]

RIE and BTE colours

Sparkling Silver · Sand · Gold · Bronze · Warm Grey · Graphite · Espresso · Deep Black

Custom colours

Anthracite · Light Beige · Beige · Medium Brown · Brown · Dark Brown

Charging options



Premium charger
Portable with powerbank
Compatible with the microRIE (CFF 8)



Standard charger
Portable
Compatible with the microRIE (CFF 8)



Desktop charger
Easy-to-use
Compatible with microRIE (CFF 8), 72 BTE (CFF 70), 76 BTE (CFF 70), 86 BTE (CFF 80)










Custom charger
Portable for ITC/ITE (CFF C)

1 Groth & Cui (2023)
2 GN Proprietary data on file

* All-day battery on one charge is based on an average use time of 20 hours, streaming 50% of the time
** Also compatible with Beltone Serene rechargeable ITC/ITE
*** Only available for ITC
**** Only available for ITE

For mild to profound hearing loss

Style	Features	Fitting range	Max output/gain°
 microRIE 62S-DRWC Rechargeable Available as CROS Transmitter	• TapControl with microRIE • All-day battery on one charge* • Receiver options: MM, LP, MP, HP and UP • Programmable push button or volume control • IP68-rated • Telecoil in 64 • CROS: CSE162S-DRWC		MM: 115/53 LP: 114/54 MP: 116/59 HP: 120/63 UP: 130/73 (Encased)
 RIE 63-DRW 312 battery RIE 64-DRW 13 battery			
 BTE 72-DWC Rechargeable  BTE 76-DWC Rechargeable Available as CROS Transmitter**	• All-day battery on one charge* • Programmable push button and/or volume control • Optional metal hook for 86-DWC • IP68-rated • Telecoil in BTE 72, 76 and 86 • DAI in BTE 86 • Exchangeable mic filters • CROS: CSE176-DWC		BTE (open): 122/51 BTE (closed): 125/58 BTE (open): 127/53 BTE (closed): 128/60
 Power BTE 86-DWC Rechargeable			Power BTE (plastic hook): 133/67 Power BTE (metal hook): 133/73
 ITC-DWC / ITE-DWC Rechargeable	• All-day battery on one charge* • Receiver options: LP***, MP, HP and UP • Shell options: ITC, Half-shell**** and Full-shell**** • Programmable push button • IP68-rated		LP: 114/40*** MP: 118/50 HP: 121/60 UP: 130/70
 ITC-DW 312 battery	• Receiver options: LP, MP, HP and UP • Shell options: ITC, Half-shell and Full-shell • Programmable push button and optional volume control • IP68-rated		LP: 114/40 MP: 117/50 HP: 120/60 UP: 130/68
 CIC non-wireless 10A battery	• Cabled fitting required • Receiver options: LP, MP and HP • Optional push button • IP68-rated		LP: 114/40 MP: 116/50 HP: 121/60

° ANSI S3.22-2014 / IEC 60118-0:2022 / 2cc coupler

Introducing **hearing clarity** and connectivity for all

Beltone **Commence**
Beltone **Boost™ SP**

Introducing Beltone Commence hearing aids for best-in-class sound quality, wearing comfort and Bluetooth® connectivity.

Best-in-class sound quality[1-4]

Best wearing comfort in tiny designs[5]

Clearest Bluetooth experience ever[6]

Colours

Shimmering miniRIE colours*

Sparkling Silver · Sand · Bronze · Espresso · Deep Black

Matte RIE colours**

Beige · Medium Blonde · Dark Brown · Sterling Grey · Black

Custom colours

Anthracite · Light Beige · Beige · Medium Brown · Brown · Dark Brown · Blue/red/ Clear shells

Charging options



Premium charger
Portable with powerbank
Compatible with the miniRIE (CFF 9)








Desktop charger
Easy-to-use
Compatible with the miniRIE (CFF 9)

1 Groth (2016)
2 Groth & Cui (2023)
3 Haastrup (2012)
4 Sjolander et al. (2019)
5 GN data on file
6 Bluetooth SIG

* Available for miniRIE 62
** Available for RIE 63, RIE 64 and Beltone Boost SP 95
*** All-day battery on one charge is based on an average use time of 20 hours, streaming 50% of the time

For mild to profound hearing loss

Style	Features	Fitting range	Max output/gain°
 miniRIE 62-DRWC Rechargeable  RIE 63-DRW 312 battery  RIE 64-DRW 13 battery	• Receiver options: LP, MP, HP and UP • All-day battery on one charge* • IP68-rated • Exchangeable mic filters • Telecoil in 64 • Programmable push button or toggle	 	LP: 113/56 MP: 115/60 HP: 120/60 UP: 129/71 (Encased)
 ITC-DW 312 battery	• Receiver options: LP, MP, HP and UP • Shell options: ITC, Half-shell and Full-shell • IP68-rated • Programmable push button and optional volume control • Bluetooth LE Audio & Auracast		LP: 113/40 MP: 115/50 HP: 120/60
 CIC non-wireless 10A battery	• Cabled fitting required • Receiver options: LP, MP and HP • IP68-rated • Optional push button		LP: 113/40 MP: 115/50 HP: 120/60
Beltone Boost SP  Super Power 95-DW 675 Battery	• Push button and toggle • Telecoil • IP68-rated • Plastic hook		Super Power BTE: 141/83

° ANSI S3.22-2014 / IEC 60118-0:2022 / 2cc coupler

Receivers, thin tubes and domes overview

	Hearing aids	Receivers and thin tubes	Domes
RIE styles	Beltone Illuminate, Beltone Envision, Beltone Commence, and Beltone Serene  **Beltone Illuminate, Beltone Envision and Beltone Serene** microRIE 62S Rechargeable **Beltone Commence** miniRIE 62 Rechargeable  RIE 63 312 battery  RIE 64 13 battery ☑ AllFit domes ☑ SureFit domes	SureFit 3 Receivers Available in 5 lengths: 0 to 4  MM* (M&RIE)  LP  MP  HP  UP (encased)	AllFit domes  Open (XS, S, M, L, XL)  Vented (XS, S, M, L, XL)  Power (S, M, L, XL)  Tulip (S, M)  Sleeve** (S, M, L)
BTE styles	Beltone Serene  **Beltone Serene** BTE 72 / BTE 76 Rechargeable ☑ AllFit domes ☑ SureFit domes	SureFit 3 thin tubes Available in 5 lengths: 0 to 4 	
Power BTE styles	Beltone Serene  **Beltone Serene** Power BTE 86*** Rechargeable		
Super Power BTE styles	Beltone Boost Max S / Beltone Boost SP  Super Power BTE 95 Metal hook Rechargeable  Super Power BTE 95 Plastic hook Rechargeable  Super Power BTE 95 Plastic hook 675 battery		

* MM is only available for Beltone Illuminate, Beltone Envision and Beltone Serene
** Recommended for M&RIE receivers
*** Metal hook option available

Our In-clinic **solutions**

First to feature NAL-NL3

Beltone Solus Max is the world's first fitting software to feature NAL-NL3[1], allowing for more personalised, data-driven fittings, alongside our proprietary fitting protocol BAFA.

Exchangeable mic filters

Reduce up to 40% of service returns with ex-changeable mic filters on all our latest RIE, Super Power BTE, and ITC styles[1]. Replacement is quick, easy, and your clients can even do it at home.

AllFit domes

Cover 97% of ear types with our lineup of AllFit domes[1], designed for the best comfort and acoustic performance[1]. Compatible with all Beltone hearing aids, including SureFit 3 and SureFit 2 receivers and thin tubes.

Aligned receivers and thin tubes

Our receivers and thin tubes are aligned to SureFit 3 across our full range of RIE and BTE hearing aids, so you do not have to work with multiple assortment boxes in your clinic.

Flexible demo devices

Our new range of demo styles lets you change technology levels easily within 10 seconds[1], so you can let your clients hear and experience the difference advanced sound features make.

1 GN Proprietary data on file (2025)

Connectivity
and accessories

Smartphone app and Remote Control

Empower your clients to take control

Download on the App Store

GET IT ON Google Play



Beltone HearMax app



Beltone Remote Control

Beltone HearMax™ app
Access to settings at all times using a mobile device

Beltone Remote Control
Control at their fingertips

Phone calls

Hands-free calls and streamed audio



Hands-free calls
Connecting hands-free to iPhone, iPad*, and Android devices compatible with Bluetooth 5.3

Direct audio streaming
From iPhone, iPad or Android devices

Multi-Mic+

Help your clients hear through the noise



Beltone Multi-Mic+

Beltone Multi-Mic+
Cut through the noise with the world's first remote mic with Auracast, now with a replaceable battery*

TV

Improve TV time for your clients



Beltone TV-Streamer+



Beltone TV Link 2


Beltone TV-Streamer+
For immediate TV speech clarity at their preferred volume. Supports Auracast broadcast audio*

Beltone TV Link 2
Enjoy crystal clear TV sound

* Compatible with Beltone Illuminate, Beltone Envision, Beltone Boost Max S, Beltone Commence and Beltone Serene (except CIC styles)


Made for **Auracast**

With Auracast broadcast audio, you can let your clients experience a new world of connectivity in public and at home.



Auracast Assistant

With the world's first Auracast Assistant for hearing aids*, every phone that is compatible with the latest version of the Beltone HearMax app can connect Auracast-compatible hearing aids to streams at train stations, airports, museums, places of worship, and more.



Multi-Mic+

Give your clients extra support in situations where they struggle to hear – over distance or in environments with unfavourable acoustics. Offering a -16dB SNR, which is a 7.7dB SNR improvement to hearing in noise compared to our leading 4-mic beamformer technology[1].

CES Innovation Awards 2025 Honoree







TV-Streamer+

The TV-Streamer+ gives an improvement of 56% in speech recognition when listening to the TV compared to listening via the TV speakers[1]. The TV-Streamer+ connects easily to help your clients listen at the volume they prefer and can transmit an Auracast stream.



1 Jespersen et al. (2024)
* GN Proprietary data on file

Feature overview

Features	Benefits	Beltone Illuminate 17	Beltone Illuminate 9	Beltone Illuminate 6	Beltone Illuminate 4	Beltone Envision 17	Beltone Envision 9	Beltone Envision 6	Beltone Envision 4	Beltone Boost Max S 17	Beltone Boost Max S 9	Beltone Boost Max S 6	Beltone Boost Max S 4	Beltone Commence & Beltone Boost SP° 4	Beltone Commence & Beltone Boost SP° 3	Beltone Commence & Beltone Boost SP° 2
Speech Intelligibility & Awareness		★★★★★	★★★★	★★★	★★	★★★★	★★★★	★★★	★★	★★★★	★★★★	★★★	★★	★★	★	★
All-Around Directionality	Better hearing with less effort in all daily life situations without the need to manually adjust the devices.	True Directionality* AI	Crosslink Directionality 4	Crosslink Directionality 3	Crosslink Directionality 2	Crosslink Directionality 4	Crosslink Directionality 3	Crosslink Directionality 3	Crosslink Directionality 2	Crosslink Directionality 4	Crosslink Directionality 3	Crosslink Directionality 3	Crosslink Directionality 2	Crosslink Directionality 2	Synchronized Speech Spotter Basic	
Speech Intelligibility in Noise	Users can control their sound experience to focus on what they want to listen to in high noise environments.	True Focus*/ Clear Focus AI		Ultra Focus II	Ultra Focus	True Focus*/Clear Focus AI	Ultra Focus II	Ultra Focus	Smart Beam Steering	Clear Focus	Ultra Focus II	Ultra Focus	Smart Beam Steering	Smart Beam Steering		Adaptive Directionality
Spatial Awareness	Vivid spatial awareness with a natural sound quality.	M&RIE/Personal Sound ID				M&RIE/Personal Sound ID				Personal Sound ID				Personal Sound ID		
Comfort		★★★★★	★★★★	★★★	★★	★★★★★	★★★★	★★★	★★	★★★★★	★★★★	★★★	★★	★★	★	★
DNN Noise Reduction *(True Noise Reduction)*	A reduction in known noises which a DNN-based system has been trained to remove for a more comfortable listening experience.	Mild – Strong* AI	On/Off* AI			Mild – Strong* AI										
Noise Reduction *(Sound Cleaner Pro)*	Improves listening comfort without compromising speech understanding.	Per environment (-10dB)	Strong (-10dB)	Considerable (-8dB)	Moderate (-6dB)	Per environment (-10dB)	Strong (-10dB)	Considerable (-8dB)	Moderate (-6dB)	Per environment (-10dB)	Strong (-10dB)	Considerable (-8dB)	Moderate (-6dB)	Moderate (-6dB)	Mild (-3dB)	
Impulse Noise Reduction	Improves comfort and sound quality for sudden sharp sounds, like clanking silverware or keys.	Mild – Strong		Mild – Moderate		Mild – Strong		Mild – Moderate		Mild – Strong		Mild – Moderate		Mild – Moderate	Mild setting	
Wind Noise Reduction	A natural and comfortable sounding experience when outdoors without the intrusion of wind turbulence noise.	Mild – Strong		Mild – Moderate		Mild – Strong		Mild – Moderate		Mild – Strong		Mild – Moderate		Mild – Moderate	Mild setting	
Silencer	Remain undisturbed by soft everyday noises such as a humming refrigerator.	Mild – Strong		Mild – Moderate		Mild – Strong		Mild – Moderate		Mild – Strong		Mild – Moderate		Mild – Moderate	Mild setting	
Acclimatization	Provides a smooth onboarding with a rising device gain over time.	Synchronized Satisfy				Synchronized Satisfy				Synchronized Satisfy				Synchronized Satisfy		
Automatic Volume Adjustments	Optimal audibility and listening comfort even when moving through rapidly changing sound environments, and without the need to make manual adjustments.	Smart Gain Pro	Smart Gain			Smart Gain Pro	Smart Gain			Smart Gain Pro	Smart Gain			Smart Gain		
Sound Processing		★★★★★	★★★★	★★★	★★	★★★★★	★★★★	★★★	★★	★★★★★	★★★★	★★★	★★	★★	★	★
Environmental Classification & Feature Engagement System	Senses where you are and provides the right support at the right time automatically.	AcoustIQ AI				AcoustIQ AI				AcoustIQ AI				AcoustIQ AI		
Warp Bands	Experience a natural and clear sound picture.	17	14	12		17	14	12		17	14	12		12	10	
Feedback Management	User is not bothered by feedback while still maintaining audibility.	Feedback Eraser 2 AI				Feedback Eraser 2 AI				Feedback Eraser 2 AI				Feedback Eraser 2 AI		
Low Frequency Boost	Boosts the listening experience with a louder, fuller sound quality when needed or preferred.	Mild – Strong				Mild-Strong				Mild-Strong				Mild-Strong		
Sound Shifter	Improves audibility of speech cues that would otherwise be lost.	Very Mild – Strong				Very Mild – Strong				Very Mild – Strong				Very Mild – Strong		
Music Processing	Preserved the rich sound quality of music.	Music Mode				Music Mode				Music Mode				Music Mode		
CROS/BICROS	A high sound quality and open sound picture for those with single-sided deafness.	Available as rechargeable microRIE style				Available as rechargeable microRIE style										
Tinnitus Sound Therapy	Personalized therapeutic sounds to relieve bothersome tinnitus.	Tinnitus Breaker Pro				Tinnitus Breaker Pro				Tinnitus Breaker Pro				Tinnitus Breaker Pro		

Connectivity

- MFi & Direct Audio Streaming with Android
- Hands-Free Calls**
- Bluetooth LE Audio & Auracast***
- Beltone HearMax app
- Streaming Bass Boost
- Multi-Mic+ & TV-Streamer+***

In-clinic ecosystem

Fitting Experience	Beltone Solus Max 2.4	Beltone Remote Care & Beltone Remote Care Live	Flexible Demo Devices***	Personalisable Patient Experience Levels	NAL-NL3 & BAFA	4 Customizable Programs
Serviceability	All-day Battery Life****	All-Weatherproof with IP68 Rating	Exchangeable Mic Filters	AllFit Domes	SureFit 3 Receivers & Thin Tubes	

Not all features are available on all form factors. For example, wireless streaming functionality and directionality features requiring ear to ear communication between devices are not available on CIC models. The feature set will also vary during a monaural fitting where some features require bilateral synchronization. For a detailed feature overview across styles, please refer to the datasheet or Beltone Solus Max fitting software catalogue.

° Beltone Boost SP is also available as Tech Level 17 in select markets. Please refer to the datasheet for more information.


* The DNN noise reduction technology used in the True Directionality and True Focus are only available on the Beltone Illuminate microRIE and Beltone Envision microRIE.

** Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation), iPad (10th generation)or later, with software updates iOS 15.3 and iPadOS 15.3 or later.

*** Not available for Beltone Boost SP and CIC models.

**** All-day battery on one charge is based on an average use time of 20 hours, streaming 50% of the time.

Nomenclature



Beltone product family

ILU*	Illuminate
ENV*	Envision
BMS	Boost Max S
BSP	Boost SP
COM	Commence
SER	Serene
CSE	Serene CROS transmitter

Technology level

17	Premium level technology
9	Advanced level technology
6	Select level technology
4	Essential plus level technology
3	Essential level technology
2	Entry level technology

Form factors

Receiver-in-Ear

62S	microRIE (Lithium-ion)
62	miniRIE (Lithium-ion)
63	Standard RIE (312)
64	Standard RIE (13)

Customs

CIC	Completely-in-Canal
ITC	In-the-Canal
ITE	In-the-Ear

Behind-the-ear

72	BTE
76	Standard BTE
86	Power BTE
86	High Power BTE
95	Super Power BTE

Additional

D	Dual microphones
R	Receiver-in-Ear
W	Wireless
C	Rechargeable
T	Telecoil
H	High Power BTE (with metal hook)

Receivers

LP	Low Power
MP	Medium Power
HP	High Power
UP	Ultra Power (encased)
MM	M&RIE (Microphone and receiver-in-ear)

Shells

FS	Full-shell
HS	Half-shell

* Available as CROS transmitter

Ready to experience Beltone?
Reach out to your sales rep for a demo or visit **gnhearing.com** to explore the portfolio and find out more.
We'll be right beside you every step of the way.

Worldwide Headquarters
GN Hearing A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: + 45 45 75 11 11
gnhearing.com

Australia
GN Hearing Australia Pty Ltd
Suite 1001, Level 1
1 Julius Avenue
North Ryde, NSW 2113
Australia
Tel: (free) 1800 658 955
gnhearing.com

M103016GB-26.05-Rev.A

© 2026 GN Hearing A/S. All rights reserved. Beltone is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone and iPad and are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast word mark and logos are trademarks owned by Bluetooth SIG, Inc. Dolby and the the double-D symbol are registered trademarks of Dolby Laboratories.