Beltone **Illuminate™**



CIC

Model*	ILU17CIC	ILU9CIC	ILU6CIC	ILU4CIC
Device Configurations				
Battery	10A			
Control options	Push button			
Receiver options	LP, MP & HP			
IP Classification	IP68			
Audiological Features				
Channels	17	14	12	12
Smart Gain Pro	●	–	–	–
Smart Gain	–	●	●	●
Sound Cleaner Pro	5 settings	4 settings	3 settings	2 settings
Silencer	3 settings	3 settings	2 settings	2 settings
Impulse Noise Reduction	3 settings	3 settings	2 settings	2 settings
Extended Bandwidth	●	●	●	●
Sound Shifter	●	●	●	●
Feedback Eraser 2 (w. AFX Music Mode)	●	●	●	●
SATISFY	●	●	●	●
Tinnitus Breaker Pro	●	●	●	●
Fitting Features				
Fitting Software Beltone Solus Max™ 2.4 or higher	●	●	●	●
Fully Flexible Programs	4	4	4	4
Safeguard Feedback Control	●	●	●	●
Datalogging	●	●	●	●

* Availability of models may vary by country.


© 2026 GN Hearing A/S. All rights reserved. Beltone is a trademark of GN Hearing A/S.

Beltone™


Technical Specifications

		LP	MP	HP	
		ANSI S3.22-2014 IEC 60118-0:2022 2cc coupler	ANSI S3.22-2014 IEC 60118-0:2022 2cc coupler	ANSI S3.22-2014 IEC 60118-0:2022 2cc coupler	
Reference test gain (60 dB SPL input)	HFA	31	36	41	dB
Full-on gain (50 dB SPL input)	Max. HFA	40 37	50 45	60 53	dB
Maximum output (90 dB SPL input)	Max. HFA	113 108	115 113	120 118	dB SPL
Total harmonic distortion	500 Hz 800 Hz 1600 Hz 3200 Hz	1.2 1.5 1.1 0.2	0.9 0.9 0.5 0.2	0.8 1.3 0.1 0.1	%
Equivalent input noise		22	22	22	dB SPL
Frequency range		≤200 - ≥8000	≤200 - ≥8000	≤200 - 7310	Hz
Current Drain (Quiescent / Operating) (Battery type 10)		0.53 / 0.73	0.55 / 0.84	0.53 / 0.76	mA

Maximum Output (OSPL 90)




Full-On and Reference Test Gain



LP
MP
HP

Fitting range



Patents pending.

All specifications are subject to change without notice.

Worldwide Headquarters
GN Hearing A/S
Lautrupbjerg 7
DK-2750 Ballerup, Denmark
Tel.: +45 45 75 11 11
beltone.com

CVR no. 55082715

Australia
GN Hearing Australia Pty Ltd
Suite 1001, Level 1
1 Julius Avenue
North Ryde, NSW 2113
Tel.: (free) 1800 658 955
www.gnhearing.com/en-au/

New Zealand
GN Hearing New Zealand Limited
4 Fred Thomas Drive
Takapuna, Auckland, 0622
Tel.: (free) 0800 900 126
www.gnhearing.com/en-nz/

4030790000GB-26.08-Rev.C