

ReSound CROS Transmitter

A solution for single-sided deafness




Improve hearing clarity all-around

Single-sided deafness is the loss of all or almost all usable hearing in one ear – making it difficult to hear speech from the impaired side, localize where sound is coming from and understand speech in noise.

The ReSound CROS and BiCROS systems can greatly enhance your daily life with a fuller sound picture and improved hearing clarity all around.

Best of all, our Organic Hearing™ philosophy makes sure the experience sounds natural, feels natural and connects you naturally to the world.

How it works



When fitted with a ReSound CROS/BiCROS solution, two ear-worn devices are required:

1. **CROS Transmitter**
(worn on the poorer ear)

   Collects sound from the poorer ear and sends the signal to the hearing aid on the better ear.

2. **Traditional hearing aid**
(worn on the better ear)

   Receives the CROS transmission, processes sound and presents the sound in the better ear. If the better hearing ear also has a mild to profound hearing loss, then a BiCROS solution will provide additional amplification for the hearing loss.

Styles to meet your needs



microRIE



BTE

ReSound CROS solutions are available in matching RIE and BTE rechargeable styles to our latest hearing aid solutions.

CROS Transmitter features

- Available in the industry's smallest rechargeable microRIE CROS and a rechargeable BTE CROS
- Provides all-day* battery life on one charge with typical use time of 16 hours
- The microRIE CROS is compatible with our premium, standard or desktop** chargers; the BTE CROS is compatible with our desktop charger
- Available in a variety of color options
- Easily adjust volume in the ReSound Smart 3D™ app

*All-day battery life on one charge is based on typical use time of 16 hours (both transmitter and hearing aid). Battery life can vary by use.

**Desktop charger is only available in the US for RIE styles.

Experience everyday benefits

Increased awareness of surroundings

Be more aware of sounds from your weaker hearing ear and enhance your overall sound perception and environmental awareness.

Improved confidence in conversation

Experience reduced feelings of isolation or frustration and improved confidence with communication in social settings.

Reduced auditory fatigue

Experience reduced listening effort with a CROS/BiCROS system, alleviating fatigue caused by straining to hear.





Learn how to wear your hearing aids correctly, personalize your sound experience, connect to Auracast broadcasts, and request assistance from your hearing care professional with the ReSound Smart 3D app.

Download on the App Store

GET IT ON Google Play

ReSound hearing aids are powered by Organic Hearing™, our unique approach to creating hearing solutions that sound natural, feel natural and connect you naturally to the world.

Explore more at resound.com

Manufacturer according to FDA:

GN ReSound North America
5005 Dean Lakes Blvd.
Shakopee, MN 55379
USA
1-888-735-4327
resound.com

Manufacturer according to Health Canada:

ReSound Canada
2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-888-737-6863
resound.com

© 2026 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone and iPad are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast word mark and logos are trademarks owned by Bluetooth SIG, Inc.

MK605954 Rev C 2026.06