ReSound GN

ReSound OMNIA • Made for hearing in noise

Try it. Hear it. Enjoy it.

Help more patients hear their best in noise with ReSound OMNIA™ – now available in a full line-up, including a new Mini RIE and rechargeable customs.



GN Making Life Sound Better

Exceptional speech understanding in noise

ReSound OMNIA offers the only beamforming in the industry to improve hearing in noise without feeling cut off. It's also the only hearing solution to simultaneously combine two benefits – narrow beamforming directionality for hearing in noise, plus omni directional listening – for a natural listening experience.[1]

The result? An incredible 150% improvement in speech understanding in noise.[2] Now you can help more patients hear their best in noise with a full line-up of popular styles.

[1] Jespersen et al (2022), Woodward et al (2022) & Braendgaard M. (2020)
[2] 4.36 dB SNR improvement compared to legacy products. Jespersen et al (2022)

Powered by the industry's only Organic Hearing

Organic Hearing™ is our unique approach to creating hearing solutions that sound natural, feel natural and connect your patients naturally to the world.

Sounds natural

Feels natural

Connects you naturally

Designed so patients can enjoy the natural sound of their own voice

One feature of ReSound OMNIA is preserving the natural sound of a user's own voice.[3] Your patients can join in conversations and talk without sounding loud or artificial to themselves. No additional features required.[3]

[3] Groth et al (2022)





Streaming without interruptions

The antenna system enables your patients to stream calls, music and more with robust connectivity and performance for a high-quality sound experience.

All-day power

Hearing confidence and worry-free use every day with battery power for a full day on one charge – even when streaming 100% of the time.*

* All-day battery on one charge is based on average use time of 12 hours when streaming 100% of the time across all styles. Battery life can vary by use.

Our best wearing comfort

Your patients can forget they are wearing hearing aids, whichever style they choose. Available in beautiful color options to match their personal style.

Each all-weatherproof design is IP68 rated with protective nanocoating on ALL components for indoor and outdoor use.



CIC
10A battery

Half shell/ITC
Rechargeable

Full shell
Rechargeable

Mini RIE
Rechargeable

RIE
Rechargeable,
13, 312 battery

BTE/Power BTE
Rechargeable

For mild to profound hearing loss

Anthracite
Light Beige
Beige
Medium Brown
Brown
Dark Brown

Custom models

Sparkling Silver
Champagne
Gold
Bronze
Warm Grey
Graphite
Espresso
Deep Black

RIE and BTE models

Colors

Check My Fit – now for ALL styles

Our insights show that wearing hearing aids correctly gives an average of 11 dB more feedback-free gain that is often lost with incorrect wear.[4] With Check My Fit, your patients can see if they are wearing their hearing aids correctly and hear their best from the start.

Find Check My Fit in the ReSound Smart 3D™ app.

[4] Yang et al (2022)


Download on the App Store

GET IT ON Google Play



Feature overview

FEATURES		TECH LEVEL 9	TECH LEVEL 7	TECH LEVEL 5
WARP bands		17	14	12
Environmental classifier		•	•	•
Environment recognition	Environmental Optimizer II	•		
	Environmental Optimizer I		•	•
Directionality	360 All-Around	STD/M&RIE		
	All Access Directionality		STD / M&RIE	
	M&RIE	•	•	•
	Spatial Sense	STD	STD	STD
	Front Focus	•		
	Ultra Focus		•	
	Binaural Directionality III			•
	Synchronized Soft Switching	•	•	•
	Autoscope	•		
	Multiscope		•	
	Adaptive Directionality			•
	Fixed Directionality	•	•	•
	Adjustable Directional Mix	•		
Noise reduction	Noise Tracker II	5 settings	3 settings	2 settings
	Personalized noise reduction	•		
Expansion		3 settings	2 settings	On / Off
Impulse Noise Reduction		3 settings	3 settings	On / Off
Wind Guard		3 settings	2 settings	On / Off
Adjustable time constants		•	•	
Multiple amplification strategies		•	•	•
Feedback management	Type	DFS Ultra III	DFS Ultra III	DFS Ultra III
	Level	3 settings	3 settings	3 settings
	Music Mode	•	•	•
	Auto DFS	•	•	•
Synchronized Acceptance Manager		•	•	•
Tinnitus Sound Generator		•	•	•
Low Frequency Boost		3 settings	2 settings	On / Off
Sound Shaper		•	•	•
e2e communication		•	•	•
Programs		4	4	4
Data Logging		•	•	•
In Situ Audiometry (not available in customs)		•	•	•
Automatic Receiver Detection		•	•	•



Manufacturer according to FDA:
GN ReSound North America
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-248-4327
pro.resound.com

ReSound Government Services
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-392-9932
gs.resound.com

Manufacturer according to Health Canada:
ReSound Canada
2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-888-737-6863
pro.resound.com

© 2023 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC.

MK605603 Rev C 2023.01