

Enhanced fitting experience

Your fitting experience matters. We're introducing a suite of improvements designed to make your everyday fittings easier and faster. With updates to both software and hardware, you can spend your time in the clinic how you want.

Exchangeable mic-filters

Aligned hearing aid basics

We're expanding the options for servicing our hearing aids in-clinic so your clients can continue to enjoy their hearing day after day.



We're aligning domes and receivers to SureFit 3 across our full new range of hearing aids[1].

To enable you to extend your service offering in-clinic, we also introduce new exchangeable mic filters on all our new Receiver-in-Ear (RIE) and Super Power Behind-the-Ear (SPBTE) styles.


Flexible demo devices

Let your clients experience the difference

Show – don't just tell.
Introducing our new range of demo styles that lets you change technology levels easily with 2 clicks[2], so you can let your clients hear and experience the difference advanced sound features make[3].

The flexible demo device allows your clients to explore various tech levels without the need to replace their physical hearing aids. This not only enhances your clients' experience but also reduces your workload - there is no need to search for new hearing aids or manage a large inventory.





- Choose preferred tech-level
- Apply to hearing aids
- Let your clients hear the rest for themselves

1 M&RIE is only available for ReSound Vivia™
2 GN Proprietary data on file
3 Flexible demo devices available in Receiver-in-Ear (RIE) and Super Power Behind-the-Ear (SPBTE) styles

Updated ReSound Smart Fit



Enjoy a smarter and more intuitive software

We are on the journey to make ReSound Smart Fit™ your software, so it works around your preferred way of working. Your feedback has been integral to ReSound Smart Fit 2.2 with improvements to the interface, faster access to onboarding reports, and even more benefits for you to explore.

The journey to faster fittings

Since we launched ReSound Smart Fit 2.0, we've seen an average of 1 minute and 19 seconds saved in fittings![1] ReSound Smart Fit is not just beautiful – it's efficient, leaving you more time to spend on client care.

The improvement to our fitting software is a journey, and with the update to ReSound Smart Fit 2.2 we've taken the next steps in creating an even more intuitive and easy fitting experience for you.

With ReSound Smart Fit 2.2 you get:

- Improved UI for an effortless flow
- Detect and fit Bluetooth®-connected hearing aids
- Accessible and improved onboarding with patient report information in Noah
- Easy phone compatibility check through QR-code scanning directly in the software
- Update both hearing aids simultaneously in less time than before, with improved and streamlined firmware updates[1]
- Our most comprehensive setup for standard verification and directionality testing, with a process that's more seamless and faster than ever[1]



Worldwide Headquarters
GN Hearing A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: +45 4575 1111
gnhearing.com

CVR no. 55082715

1 GN Proprietary data on file

© 2025 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S.
The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.

M201716GB-25.06-Rev.A