Hear through noise without effort

Beltone **Achieve**™

Experience more of life with the outstanding solution for speech understanding in noise.*



Beltone™

* Jespersen et al (2022)

An outstanding solution and **a full family**

Offer your clients outstanding speech understanding in noise[1], without effort, with the full family of Beltone Achieve hearing aids, now including a miniRIE. You will experience positive impacts on your practice and customer loyalty – greater client satisfaction, improved retention and fewer revisits and returns.











Increase client satisfaction and wear time

Hear through the noise

With Beltone Achieve you can tackle your clients' biggest hearing challenge – speech understanding in noise – and give a great fit for this from the start. Offer our best speech-understanding-in-noise technology[1], use our fitting software optimisations, and look forward to successful first fits for greater retention and fewer follow-ups.

Greater uninterrupted streaming

Let your clients enjoy clearer sound[3] and seamlessly stream calls, music, and more with double the improvement and an extra-robust antenna[3]. This can motivate them to wear their hearing aids for longer, to increase satisfaction and reduce the need for service and support.

Made for iPhone | iPad | iPod

Works with android

Support each client's own voice

Natural own voice perception lets your clients join in conversations without the occlusion effect, so their voices don't sound loud or artificial to them through the hearing aids. This helps ensure a greater first fit acceptance, to increase retention.[2]

Confidence in any moment, rain or shine

Give a worry-free experience with all-day battery power on 1 charge – even when streaming 100% of the time.* The durable Beltone Achieve designs are all-weatherproof and IP68 rated with nanocoating on every single part.[4]

Hearing comfort throughout the day

All Beltone Achieve styles are designed for comfort, so your clients can forget they're wearing hearing aids. Our dedication to care and comfortable hearing leads to greater satisfaction and fewer returns.



CIC
10A battery

ITE/ITC
Rechargeable

MiniRIE
Rechargeable RIE

RIE
Rechargeable,
13, 312 battery

BTE
Rechargeable BTE
and Power BTE

For mild to profound hearing loss

Custom Colors

- Anthracite
- Beige Light
- Beige
- Brown Medium
- Brown
- Brown Dark

Behind the Ear colors

- Sparkling Silver
- Champagne
- Gold
- Bronze
- Warm Grey
- Graphite
- Espresso
- Deep Black

1. Jespersen et al (2022)
2. Data on file (2022)
3. Compared to legacy technology. Data on file (2022)
4. Technical data on file (2022)

* All-day battery on 1 charge is based on average use time of 12 hours when streaming 100% of the time across all styles. Battery life can vary by use.

Feature overview

		Beltone Achieve			
		IMP17	IMP9	IMP7	IMP4
FEATURE	Wide Dynamic Range Compression	17	14	12	12
	Environmental Classifier	●	●	●	●
SPEECH IN NOISE	Ultra Focus 2	●	–	–	–
	Ultra Focus	–	●	–	–
	CrossLink Directionality 4	●	–	–	–
	CrossLink Directionality 3	–	●	–	–
	CrossLink Directionality 2	–	●	●	–
	CrossLink Directionality	–	–	–	●
	Personal Sound ID	●	●	●	●
	Spatial Directionality	–	–	–	–
	Synchronized Speech Spotter Pro	●	●	–	–
	Synchronized Speech Spotter Basic	–	–	●	●
	Speech Spotter Basic	–	–	–	–
	Smart Beam Steering	●	●	–	–
	Adaptive Directionality	–	–	●	●
	Fixed Hypercardioid	●	●	●	●
	Band-split Directionality	●	●	●	●
	Adjustable Mixing Point Frequency	●	–	–	–
HEARING COMFORT	Environmental Programs	4	4	4	4
	Environmental Gain Steering	Smart Gain Pro	Smart Gain	Smart Gain	Smart Gain
	Silencer (expansion)	3 settings	2 settings	On/Off	On/Off
	Sound Cleaner (noise reduction)	5 settings	3 settings	2 settings	On/Off
	Impulse Noise Reduction	3 settings	2 settings	On/Off	On/Off
	Wind Noise Reduction	3 settings	2 settings	On/Off	On/Off
	Low Frequency Boost	3 settings	2 settings	On/Off	On/Off
	Safeguard Feedback Control	●	●	●	●
	Feedback Eraser 2	●	●	●	●
	Feedback Eraser	–	–	–	–
	AFX Music Mode	●	●	●	●
	Amplification Strategy	●	●	●	●
	Synchronized SATISFY (adaptation management)	●	●	●	●
ADDITIONAL FEATURES	Sound Shifter (frequency lowering)	●	●	●	●
	Tinnitus Breaker Pro	●	●	●	●
	Auto-Phone	●	●	●	●
	Delayed Activation	●	●	●	●
	Datalogging (Satisfaction Journal)	●	●	●	●
FUNCTIONAL FEATURES	Hands-free phone and FaceTime calls for iPhone and iPad*	●	●	●	●
	e2e communication	●	●	●	●
	Asymmetric Phone Handling	●	●	●	●
	Wireless accessories	●	●	●	●
	Direct audio streaming from iPhone, iPad or iPod (MFi)	●	●	●	●
	Direct audio streaming from Android™ devices	●	●	●	●
	Beltone HearMax™ app	●	●	●	●

*Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later.

Ready to experience Beltone Achieve?

Reach out to your sales rep for a demo or visit beltone.com to find out more.
We'll be right beside you every step of the way.

Worldwide Headquarters
Beltone A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: + 45 45 75 11 11
beltone.com

Australia
GN Hearing Australia Pty Ltd
Gate C, 19-25 Khartoum Rd
Macquarie Technology Park
Macquarie Park NSW 2113
Tel.: (free) 1800 658 955
beltone.com/en-au

New Zealand
GN Hearing New Zealand Limited
Ground Floor, North Entrance
4 Fred Thomas Drive
Takapuna, Auckland, 0622
Tel.: (free) 0800 900 126
beltone.com/en-au

Beltone™

© 2023 GN Hearing A/S. All rights reserved. Beltone is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone, iPad, iPod touch and FaceTime are trademarks of Apple Inc., registered in the U.S. and other countries. Android is a trademark of Google LLC.

M102909GB-22.11-Rev.A